

One Step Multi-Drug Urine Test, 12 A



Drugs are substances that affect the structure or function of the body and are taken with the intention of diagnosis, treatment, or prevention of disease.

Drug abuse (or substance abuse) refers to the use of drugs without medical supervision, usually for their mood and behavior altering effects.

The use of psychoactive drugs is prohibited or controlled by governing authorities; due to its devastating health effects, higher morbidity and mortality risks, potential for dependency and addiction, negative social impacts; Its worth noting that the rate of overdose fatalities increased during the COVID-19 pandemic.

INTENDED USE

Turklab's One Step Multi-Drug Urine Test is a rapid urine screening immunochromatographic assay, that can be performed without the use of an instrument, for the qualitative detection of multiple drugs and drug metabolites in urine at the following cut-off concentrations:

PARAMETERS	CUT-OFF LEVELS
Amphetamine (AMP)	500 ng/mL
Barbiturates (BAR)	300 ng/mL
Benzodiazepines (BZO)	300 ng/mL
Cocaine (COC)	150 ng/mL
Marijuana (THC)	50 ng/mL
Methadone (MTD)	300 ng/mL
Methamphetamine (MET)	500 ng/mL
Methylenedioxymethamphetamine (MDMA)	500 ng/mL
Opiate (OPI)	2000 ng/mL
Tricyclic Antidepressants (TCA)	1000 ng/mL
Buprenorphine (BUP)	10 ng/mL
Synthetic Marijuana (K2)	50 or 25 ng/mL



PRODUCT	TEST TYPE	SAMPLE TYPE	KIT SIZE	SENSITIVITY	READING TIME	KIT COMPONENTS
ONE STEP MULTI-DRUG URINE TEST	Panel	Urine	25 Test Panel	Variable	5 Min.	Test Panel, Instructions For Use